

## KIT DE RÉTROFIT AVEC DOUCHETTE POUR DOUCHE DE SÉCURITÉ ET LAVEUR D'YEUX HAWS

### RETROFIT DRENCH HOSE KIT FOR HAWS SAFETY SHOWERS AND EYEWASHES




## CARACTÉRISTIQUES / SPECIFICATIONS

Le Haws EU-8901RFK-AXION est un kit de retrofit avec douchette permettant d'ajouter un dispositif de rinçage manuel à une douche de sécurité ou à un laveur d'yeux existant. Il permet d'effectuer un rinçage ciblé du corps, des yeux ou du visage tout en conservant l'installation déjà en place.

La douchette intègre la technologie AXION® de Haws, conçue pour assurer une distribution uniforme et douce de l'eau. Elle permet de diriger le rinçage précisément vers la zone concernée en cas d'exposition à des produits chimiques ou à des substances dangereuses.

Conçu pour l'adaptation d'équipements existants, le kit EU-8901RFK-AXION comprend les différents éléments nécessaires à son installation, notamment la douchette AXION®, les colliers de serrage, les raccords en T, les raccords et le matériel de fixation. Il permet ainsi de compléter une installation de sécurité sans avoir à remplacer l'ensemble de l'équipement.

Ce kit convient notamment aux sites industriels, laboratoires, ateliers, plateformes logistiques et autres environnements professionnels utilisant des douches de sécurité ou des laveurs d'yeux Haws et souhaitant disposer d'un dispositif de rinçage manuel supplémentaire.



Free access to the eye/facewash & body spray EU-8901RFK-AXION must be guaranteed at all times.  
The operator is responsible for the observance of this requirement.



# INSTALLATION, FONCTIONNEMENT ET MAINTENANCE / INSTALLATION, OPERATION AND MAINTENANCE INSTRUCTIONS

## Shower location :

The eye/facewash & body spray EU-8901RFK-AXION must be installed in the vicinity of potential danger zones. It must be highly visible and easily accessible.

## Connection conduit :

You are recommended to use a 1/2" G conduit, capable of withstanding pressures of 2 - 6BAR (30 - 90PSI). The use of an inlet filter designed to prevent the entry of suspended particles and other impurities is also recommended.

## Water quality :

Eye/facewash & body spray require a supply of drinking water conforming to European standards or their equivalents in other countries.

## Water temperature :

The required water temperature should be determined by the user. The accepted ideal range is 15 - 35°C / 60 - 95°F (but please ensure conformity with local regulations).

## Volume flow :

A minimum flow rate of 6 litres / 1,6 US gallons per minute must be maintained. The eye/facewash & body spray must be capable of delivering this amount of water for at least 15 minutes.

## Connection to plumbing :

The inlet consists of a 1 1/4" G male-threaded conduit.

## Regular inspection :

The eye/facewash & body spray 8901RFK Axion should be inspected at weekly intervals. The inspection tag should be dated and signed after each check.

| TROUBLESHOOTING               |                                                                                                                                                                                                                                  |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| PROBLEM                       | REPAIR CHECKLIST                                                                                                                                                                                                                 |
| 1. No flow                    | Check the main shutoff valve.                                                                                                                                                                                                    |
| 2. Insufficient flow of water | 1. Check the line pressure. This should be at least 2 BAR / 30 PSI (dynamic).<br>2. Blocked flow, possibly due to the conduit not being correctly flushed out. Dismantle the spray element, and clean the rubber flow-regulator. |